

NEO Shield- Series™

Creating New Dimensions for Neonatal & Child Care



PRODUCT CATALOGUE



www.neokraftmedical.com

 **NEOKRAFT™**
ISO 13485 Certified Company

Bubble CPAP System



NEOPAP - 505



Air / Oxygen Mixer with high pressure protection

Servo Humidifier with invasive & non-invasive modes

Low noise medical grade Air Compressor

Infant heated wired breathing circuit to maintain air way temperature

Bubble Generator

Mounted on anti skid casters with locking facility

Model NEOPAP-505

Range of Respiratory Humidifiers



Servo Humidifier



Non-Servo Humidifier with Display



Non-Servo Humidifier

* Technical specifications may be altered by Neokraft without prior notice, accessories shown in the catalogue are not part of the standard equipment

Manufactured by:



NEOKRAFT MEDICAL PRIVATE LIMITED

18, 1st Cross, 11th Main, Tata Nagar, Sahakara Nagar [P],
BANGALORE - 560 092, INDIA.

Phone : (080) 6500 0249

Mobile : +91 9686 78 5555

Email : info@neokraftmedical.com

Website : www.neokraftmedical.com

Channel Partner